

ZNF557 Rabbit Polyclonal Antibody



CAB18490

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

49kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

79230

Uniprot

Q8N988

Synonyms:

ZNF557

Immunogen:

Recombinant protein of human ZNF557.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZNF557 Rabbit pAb (CAB18490)