

SLC52A2 Rabbit Polyclonal Antibody



CAB18491

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79581

Uniprot

Q9HAB3

Synonyms:

BVVLS2; D15Ert747e; GPCR41; GPR172A; PAR1; RFT3; RFVT2;
hRFT3; SLC52A2

Immunogen:

Recombinant protein of human SLC52A2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
