

# ARHGAP10 Rabbit Polyclonal Antibody



CAB18492

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

79658

#### Uniprot

A1A4S6

#### Synonyms:

GRAF2; PS-GAP; PSGAP; ARHGAP10

#### Immunogen:

Recombinant protein of human ARHGAP10.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---