

Product Information

Product SKU:	CAB18493	Gene ID:	79811	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	117kDa	Conjugate:	-
Calculated MW:	117kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables RNA binding activity. Predicted to be involved in regulation of mRNA processing and regulation of transcription by RNA polymerase II. Located in nuclear body.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	Met; SLTM
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human SLTM (NP_079031.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.