

CAB18499

Product Information

Product SKU:	CAB18499	Gene ID:	80325	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	-
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein with an ankyrin repeat region and two BTB/POZ domains, which are thought to be involved in protein-protein interactions. Expression of this gene is activated by the phosphatase and tensin homolog, a tumor suppressor. Alternate splicing results in three transcript variants.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	BPOZ; BTB3; BTBD21; EF1ABP; PP2259; ABTB1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ABTB1 (NP_742024.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.