

ABTB1 Rabbit Polyclonal Antibody



CAB18499

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

55kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

80325

Uniprot

Q969K4

Synonyms:

BPOZ; BTB3; BTBD21; EF1ABP; PP2259; ABTB1

Immunogen:

Recombinant protein of human ABTB1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - ABTB1 Polyclonal Antibody (CAB18499)

