

Product Information

Product SKU:	CAB18522	Gene ID:	91947	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	45kDa	Conjugate:	-
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquitin ligase-substrate adaptor activity. Acts upstream of or within positive regulation of ubiquitin-protein transferase activity. Located in endosome and plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ARRDC4
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ARRDC4 (NP_899232.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.