

CAB18528

Product Information

Product SKU:	CAB18528	Gene ID:	114781	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	-
Calculated MW:	69kDa	Isotype:	IgG

Immunogen Information

Background: This locus encodes a BTB/POZ domain-containing protein. This domain is known to be involved in protein-protein interactions. Polymorphisms at this locus have been reported to be associated with susceptibility to Restless Legs Syndrome and may also be associated with Tourette Syndrome. Alternatively spliced transcript variants have been described.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: dJ322I12.1; BTBD9

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 283-582 of human BTBD9 (AAI01358.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.