

BTBD9 Rabbit Polyclonal Antibody



CAB18528

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

70kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

114781

Uniprot

Q96Q07

Synonyms:

dJ322112.1; BTBD9


Immunogen:

Recombinant protein of human BTBD9.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - BTBD9 Rabbit pAb (CAB18528)