

GTSF1 Rabbit Polyclonal Antibody



CAB18534

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

19kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

121355

Uniprot

Q8WW33

Synonyms:

FAM112B; GTSF1

Immunogen:

Recombinant protein of human GTSF1.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - GTSF1 Rabbit pAb (CAB18534)