

FAM84B Rabbit Polyclonal Antibody



CAB18539

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

34kDa

Applications:

WB

Reactivity:

Human

Protein Background

Immunogen information

Gene ID:

157638

Uniprot

Q96KN1

Synonyms:

NSE2; LRATD2; BCMP101; FAM84B

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 151-310 of human FAM84B (NP_777571.1).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
