

TMEM18 Rabbit Polyclonal Antibody



CAB18541

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

23kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

129787

Uniprot

Q96B42

Synonyms:

TMEM18


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 40-100 of human TMEM18 (NP_690047.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - TMEM18 Rabbit pAb (CAB18541)