

# C14orf37 Rabbit Polyclonal Antibody



CAB18547

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

145407

#### Uniprot

Q86TY3

#### Synonyms:

UT2; c14\_5376; C14orf37

#### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 700-774 of human C14orf37 (NP\_001001872.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---