
Product Information

Product SKU:	CAB18551	Gene ID:	151636	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	84kDa	Conjugate:	-
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background: Enables several functions, including STAT family protein binding activity; ubiquitin-like protein ligase binding activity; and ubiquitin-protein transferase activity. Involved in several processes, including positive regulation of macromolecule metabolic process; positive regulation of protein localization; and protein ubiquitination. Located in cytosol; lysosome; and nucleoplasm. Part of protein-containing complex. Colocalizes with early endosome and nucleus.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: BBAP; RNF143; DTX3L

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 601-740 of human DTX3L (NP_612144.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.