

# ZNF645 Rabbit Polyclonal Antibody



CAB18553

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

49kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

158506

#### Uniprot

Q8N7E2

#### Synonyms:

CT138; HAKAIL; ZNF645

#### Immunogen:

Recombinant protein of human ZNF645.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - ZNF645 Rabbit pAb (CAB18553)

