
Product Information

Product SKU:	CAB18553	Gene ID:	158506	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	49kDa	Conjugate:	-
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the zinc finger domain-containing protein family. This family member contains both a RING-type and a C2H2-type of zinc finger domain, and it may function as an E3 ubiquitin-protein ligase. Protein localization suggests a role in human sperm production and quality control.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: CT138; HAKAIL; ZNF645

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 346-425 of human ZNF645 (NP_689790.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.