

WBP2NL Rabbit Polyclonal Antibody



CAB18555

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

32kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

164684

Uniprot

Q6ICG8

Synonyms:

GRAMD7; PAWP; WBP2NL


Immunogen:

Recombinant protein of human WBP2NL.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - WBP2NL Rabbit pAb (CAB18555)