

LEMD2 Rabbit Polyclonal Antibody



CAB18558

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

57kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

221496

Uniprot

Q8NC56

Synonyms:

CTRCT42; NET25; dJ482C21.1; LEMD2


Immunogen:

Recombinant protein of human LEMD2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - LEMD2 Rabbit pAb (CAB18558)