

Product Information

Product SKU:	CAB18558	Gene ID:	221496	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	57kDa	Conjugate:	-
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene encodes a LEM domain-containing transmembrane protein of the inner nuclear membrane. The protein is involved in nuclear structure organization and plays a role in cell signaling and differentiation. Mutations in this gene result in Cataract 46, juvenile-onset. Multiple transcript variants have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	LEM2; NET25; MARUPS; CTRCT42; dJ482C21.1; LEMD2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human LEMD2 (NP_851853.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.