

ZNF326 Rabbit Polyclonal Antibody



CAB18562

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

284695

Uniprot

Q5BKZ1

Synonyms:

ZAN75; ZIRD; Zfp326; dJ871E2.1; ZNF326

Immunogen:

Recombinant protein of human ZNF326.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
