

Product Information

Product SKU:	CAB18568	Gene ID:	375337	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	191kDa	Conjugate:	-
Calculated MW:	191kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to be involved in spermatid development and spermatocyte division. Predicted to act upstream of or within apoptotic process; ncRNA transcription; and positive regulation of meiotic cell cycle phase transition. Predicted to be located in cytosol.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	C3orf77; TAZ1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 13-276 of mouse TAZ1 (NP_001186665.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.