

TAZ1 Rabbit Polyclonal Antibody



CAB18568

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

191KDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

375337

Uniprot

Q8N9V7

Synonyms:

C3orf77; TAZ1


Immunogen:

Recombinant protein of human TAZ1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - TAZ1 Rabbit pAb (CAB18568)