

# PPP2R1A/PPP2R2A/PPP2R1B Rabbit Polyclonal Antibody

## CAB18571



### Product Information

**Size:**

20uL, 50uL, 100uL, 200uL

**Observed MW:**

55kDa/66kDa

**Calculated MW:****Applications:**

WB

**Reactivity:**

Human, Mouse, Rat

### Antibody Information

**Recommended dilutions:**

WB 1:500 - 1:2000

**Source:**

Rabbit

**Isotype:**

IgG

**Purification:**

Affinity purification

### Protein Background

#### Immunogen information

**Gene ID:**

5518/ 5520/ 5519

**Uniprot**

P30153/P63151/P30154

**Synonyms:**

PPP2R1A/PPP2R2A/PPP2R1B


**Immunogen:**

Recombinant protein of human PPP2R1A/PPP2R2A/PPP2R1B.

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images



Western blot - PPP2R1A/PPP2R2A/PPP2R1B pAb (CAB18571)