

# PRELP Rabbit Polyclonal Antibody



CAB18573

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

44KDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

5549

#### Uniprot

P51888

#### Synonyms:

MST161; MSTP161; SLRR2A; PRELP

#### Immunogen:


Recombinant protein of human PRELP.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



^- PRELP Polyclonal Antibody (CAB18573)