

POU3F2 Rabbit Polyclonal Antibody



CAB18576

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

55KDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5454

Uniprot

P20265

Synonyms:

BRN2; N-Oct3; OCT7; OTF-7; OTF7; POUF3; brn-2; oct-7; POU3F2


Immunogen:

Recombinant protein of human POU3F2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of POU3F2 Polyclonal Antibody (CAB18576)