

nifD Rabbit Polyclonal Antibody



CAB18578

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Protein Background

Immunogen information

Gene ID:

1054811

Uniprot

Applications:

WB

Synonyms:

nifD

Reactivity:

Bradyrhizobium
diazoefficiens

Immunogen:

Recombinant protein of Bradyrhizobium diazoefficiens nifD.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
