

RPS21 Rabbit Polyclonal Antibody



CAB18585

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

13kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

6227

Uniprot

P63220

Synonyms:

HLDF; S21; RPS21

Immunogen:

Recombinant protein of human RPS21.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - RPS21 Rabbit pAb (CAB18585)

