

Product Information

Product SKU:	CAB18599	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	12kDa/30kDa	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:	Protein N-myristoylation is an important fatty acylation catalyzed by N-myristoyltransferases (NMTs), which are ubiquitous enzymes in eukaryotes. Specifically, attachment of a myristoyl group is vital for proteins participating in various biological functions, including signal transduction, cellular localization, and oncogenesis.
Recommended Dilution:	WB, 1:100 - 1:500
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence containing myristylated G.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.