

PI3 Kinase p85 alpha Mouse Monoclonal Antibody

CAB18640



Product Information

Size:

50uL

Observed MW:

84kDa

Calculated MW:**Applications:**

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Mouse

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5295

Uniprot

P27986

Synonyms:

PI3K; p50alpha; p55alpha; p85alpha


Immunogen:

Recombinant protein of Human PI3 Kinase p85 alpha.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PI3 Kinase p85 alpha Mouse mAb (CAB18640)