

# TM6SF2 Rabbit Polyclonal Antibody



CAB18649

---

## Product Information

### Size:

50uL, 100uL, 200uL

### Observed MW:

43kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse,

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

53345

#### Uniprot

Q9BZW4

### Synonyms:

#### Immunogen:


Recombinant protein of human TM6SF2.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - TM6SF2 Rabbit pAb (CAB18649)