

APOBEC2 Rabbit Polyclonal Antibody



CAB18653

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

26kDa

Calculated MW:

25kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10930

Uniprot

Q9Y235

Synonyms:

APOBEC2; ARCD1; ARP1; C->U-editing enzyme APOBEC-2

Immunogen:

Recombinant protein of human APOBEC2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - APOBEC2 Rabbit pAb (CAB18653)

