
Product Information

Product SKU:	CAB18657	Gene ID:	22834	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	72kDa	Conjugate:	-
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human ZNF652 (NP_055712.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.