

H2-Ab1 Rabbit Polyclonal Antibody



CAB18658

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

30kDa

Calculated MW:

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

14961

Uniprot

P06345

Synonyms:

A1845868; Abeta; H-2Ab; H2-Ab; I-Abeta; IAb; Ia-2; Ia2; Rmcs1

Immunogen:

Recombinant protein of Mouse H2-Ab1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - H2-Ab1 Rabbit pAb (CAB18658)

