

YEATS2 Rabbit Polyclonal Antibody



CAB18666

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

151kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55689

Uniprot

Q9ULM3

Synonyms:


Immunogen:

Recombinant protein of human YEATS2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - YEATS2 Rabbit pAb (CAB18666)