

PYM1 Rabbit Polyclonal Antibody



CAB18669

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

23kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

84305

Uniprot

Q9BRP8

Synonyms:

Immunogen:

Recombinant protein of human PYM1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - PYM1 Rabbit pAb (CAB18669)

