

CAB18669

Product Information

|              |          |               |        |             |                 |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB18669 | Gene ID:      | 84305  | Size:       | 20uL, 100uL     |
| Clone No:    | -        | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

|                |       |            |     |
|----------------|-------|------------|-----|
| Observed MW:   | 23kDa | Conjugate: | -   |
| Calculated MW: | 23kDa | Isotype:   | IgG |

Immunogen Information

**Background:** Enables ribosome binding activity. Involved in exon-exon junction complex disassembly; nuclear-transcribed mRNA catabolic process, nonsense-mediated decay; and positive regulation of translation. Located in cell junction; cytosol; and nuclear lumen. Part of exon-exon junction complex.

|                       |  |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200  |
| Synonyms:             | PYM; WIBG; PYM1  |
| Purification Method:  | Affinity purification  |
| Immunogen:            | Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human PYM1 (NP_115721.1). |
| Storage:              | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.                 |