

SMAD2/SMAD3 Rabbit Polyclonal Antibody



CAB18674

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

52kDa/60kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4087/4088

Uniprot

Q15796/P84022

Synonyms:


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 350 to 425 of human SMAD2/SMAD3 (NP_005893.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SMAD2/SMAD3 Rabbit pAb (CAB18674)