

LOC101736823 Rabbit Polyclonal Antibody



CAB18679

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Bombyx mori

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

101736823

Uniprot

Synonyms:

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 251-460 of Bombyx mori LOC101736823.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
