

ABCB1 Rabbit Polyclonal Antibody



CAB18684

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

150kDa/180kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

11423048

Uniprot

Synonyms:


Immunogen:

Recombinant protein of Medicago truncatula ABCB1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ABCB1 Rabbit pAb (CAB18684)