

FFAR4 Rabbit Polyclonal Antibody



CAB18689

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

42KDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

338557

Uniprot

Q5NUL3

Synonyms:

Immunogen:

Recombinant protein of human FFAR4.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - FFAR4 Rabbit pAb (CAB18689)

