

CAB18848

Product Information

Product SKU:	CAB18848	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Species independent

Additional Information

Observed MW:	25kDa(Light chain),55kDa(Heavy chain)	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: It is well know now that biotin plays important roles in a variety of critical metabolic reactions in the cell, and thus, is essential for normal human health, growth and development. Biotin is a water-soluble vitamin required by all organisms by virtue of its essential role in carboxylation reactions. Streptavidin and its homologs (together referred to as streptavidin) are widely used in molecular science owing to their highly selective and stable interaction with biotin.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Biotin
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.