

# SRD5A1 Rabbit Polyclonal Antibody



CAB1898

---

## Product Information

### Size:

50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

29kDa

### Applications:

WB IHC

### Reactivity:

Human

## Protein Background

Steroid 5-alpha-reductase (EC 1.3.99.5) catalyzes the conversion of testosterone into the more potent androgen, dihydrotestosterone (DHT). Also see SRD5A2 (MIM 607306).[supplied by OMIM, Mar 2008]

## Immunogen information

### Gene ID:

6715

### Uniprot

P18405

### Synonyms:

SRD5A1; S5AR 1; S5AR1

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant protein of human SRD5A1

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---