
Product Information

Product SKU:	CAB19001	Gene ID:	8019	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	80kDa	Isotype:	IgG

Immunogen Information

Background: This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ORFX; FSHRG2; RING3L; BRD3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human BRD3 (NP_031397.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.