

BRD3 Rabbit Polyclonal Antibody



CAB19001

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

110kDa

Calculated MW:

60kDa/79kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known.

Immunogen information

Gene ID:

8019

Uniprot

Q15059

Synonyms:

BRD3; ORFX; RING3L

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human BRD3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - BRD3 Polyclonal Antibody (CAB19001)

