

CAB19098

Product Information

Product SKU:	CAB19098	Gene ID:	5816	Size:	20uL, 100uL
Clone No:	ARC0385	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	12kDa	Conjugate:	Unconjugated
Calculated MW:	12kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a high affinity calcium ion-binding protein that is structurally and functionally similar to calmodulin and troponin C. The encoded protein is thought to be involved in muscle relaxation. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: D22S749; Parvalbumin (PVALB)

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human Parvalbumin (PVALB) (P20472).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.