

CAB19140

Product Information

Product SKU:	CAB19140	Gene ID:	4324115	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Component of the RACK1 regulatory proteins that functions in innate immunity by interacting with multiple proteins in the RAC1 immune complex. Acts as positive regulator of reactive oxygen species (ROS) production and is required for resistance against rice blast (M.grisea) infection.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RWD; GPB-LR; RACK1A; OsRACK1A
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsRACK1A.
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.