

**CAB1919**

---

**Product Information**

<b>Product SKU:</b>	CAB1919	<b>Gene ID:</b>	7001	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	22kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	22kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a member of the peroxiredoxin family of antioxidant enzymes, which reduce hydrogen peroxide and alkyl hydroperoxides. The encoded protein plays an antioxidant protective role in cells, and it may contribute to the antiviral activity of CD8(+) T-cells. The crystal structure of this protein has been resolved to 2.7 angstroms. This protein prevents hemolytic anemia from oxidative stress by stabilizing hemoglobin, thus making this gene a therapeutic target for patients with hemolytic anemia. This protein may have a proliferative effect and play a role in cancer development or progression. Related pseudogenes have been identified on chromosomes 5, 6, 10 and 13.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** PRP; TSA; PRX2; PTX1; TPX1; NKEFB; PRXII; TDPX1; NKEF-B; HEL-S-2a; 2)

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Peroxiredoxin 2 (Peroxiredoxin 2 (PRDX2)) (NP\_859428.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.