

## CAB19270

---

**Product Information**

<b>Product SKU:</b>	CAB19270	<b>Gene ID:</b>	51701	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2441	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	58kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	58kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** Enables ubiquitin protein ligase binding activity. Involved in protein stabilization and transforming growth factor beta receptor signaling pathway. Predicted to be located in cytosol and nucleoplasm. Predicted to be active in cytoplasm and nucleus.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** NLK; nemo like kinase; LAK1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 400-500 of human NLK (NP\_057315.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.