

Product Information

Product SKU:	CAB19272	Gene ID:	4144	Size:	20uL, 100uL
Clone No:	ARC2447	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

**Background:** The protein encoded by this gene catalyzes the production of S-adenosylmethionine (AdoMet) from methionine and ATP. AdoMet is the key methyl donor in cellular processes.

Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	MATA2; MATI; SAMS2; MAT2A
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAT2A (P31153).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.