

Product Information

Product SKU:	CAB19290	Gene ID:	1153	Size:	20uL, 100uL
Clone No:	ARC2472	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	19kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables mRNA 3'-UTR binding activity and small ribosomal subunit rRNA binding activity. Involved in mRNA stabilization; positive regulation of translation; and response to UV. Located in cytoplasm and nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:500 - 1:1000 IF/ICC,1:50 - 1:200

**Synonyms:** CIRP; CIRBP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CIRBP (Q14011).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.